

3.5 mm Distal Tibia With Tab  
Locking Plate




**Bonetech Medisys**  
implanting precision



**Bonetech Medisys**  
implanting precision

## 3.5 mm Distal Tibia With Tab Locking Plate



### Feature :

- Anatomically contoured, plate is twisted 20° and bent to fit the distal tibia.
- Locking holes in distal region are parallel to the joint.
- Distal tab for optional medial malleolus Screw accepts 3.5 mm locking, 2.7 mm cortex, 3.5 mm cortex.

HOLE	CODE (Left)	CODE (Right)
5	BT359.T05L	BT359.T05R
7	BT359.T07L	BT359.T07R
9	BT359.T09L	BT359.T09R
11	BT359.T11L	BT359.T11R
13	BT359.T13L	BT359.T13R
T-Titanium, S-Steel		